# **RIDASCREEN<sup>®</sup> FAST Ei / Egg Protein**

**Art. Nr. 6402**

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- Probenaufarbeitung für Wein
- RIDASCREEN®FAST Ei/Egg Protein (Art. No. R6402) used with automates GEMINI and ChemWell® 2910 Automate
- RIDASCREEN®FAST – Swabbing Allergen Swabbing method for the qualitative analysis of allergens in the production line or for laboratory equipment

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# RIDASCREEN<sup>®</sup>FAST Ei / Egg Protein

## Brief information

RIDASCREEN<sup>®</sup>FAST Ei / Egg Protein (Art. No. R6402) is a sandwich enzyme immunoassay for the quantitative analysis of whole egg (-powder) in food. All reagents required for the enzyme immunoassay - including standards - are contained in the test kit.

The test kit is sufficient for 48 determinations (including standards). A microtiter plate spectrophotometer is required for quantification.

- Sample preparation: homogenization, extraction and centrifugation
- Time requirement: sample preparation (for 10 samples) ..... approx. 20 min  
test implementation (incubation time) ..... 30 min
- Standard material: NIST Reference Material 8445 (whole egg powder)
- Limit of detection: 0.10 mg/kg (ppm) whole egg powder  
(corresponds to 0.03 mg/kg (ppm) egg white protein)
- Limit of quantification: 0.5 mg/kg (ppm) whole egg powder  
(corresponds to 0.13 mg/kg (ppm) egg white protein)
- Specificity: The specific antibodies detect the egg white proteins ovalbumin and ovomucoid.  
No cross reactivity was observed to raw or cooked meat of chicken, turkey hen, pork and beef.

Cross reactivities of the used antibodies have been determined for the pure food (e.g. corn flour). In a composed / processed food (e.g. maize bread) cross reactivities might be different. Interfering substances (e.g. polyphenols) can be detected by spike experiments.

In order to increase the quality of assessment when performing ELISA-procedures, we refer additionally to our Good ELISA Practice (GEP) – Manual in

the respective version. Such lists minimum standards concerning the framework conditions when using test-kits of R-Biopharm AG and performing ELISA-analysis. The manual can be retrieved, printed and downloaded under the website [www.r-biopharm.com/products/food-feed-analysis](http://www.r-biopharm.com/products/food-feed-analysis).

## Related products

Bioavid Lateral Flow Ei/Egg (BL608-10 & -25)

### 1. Intended use

RIDASCREEN®FAST Ei / Egg Protein is a sandwich enzyme immunoassay for the quantitative analysis of whole egg (-powder) in food like salad dressings, sausages, wines, baking-mixtures for cakes or bread and ice cream.

### 2. General

The allergen can be present as an ingredient or as a contamination in raw and cooked products. According to the **regulation (EU) No. 1169/2011** from egg and products thereof must be declared on food labels. Similar regulations exist e.g. in the USA, Canada, Australia and New Zealand.

Egg white contains 9 - 11 % of protein. Four main allergenic proteins come to 80 % of the egg white proteins. These main allergens are ovomucoid (11 %), ovalbumin (54 %), ovomucin (12 %) and lysozyme (3.5 %). The allergenic potential of the proteins present in egg yolk is only moderate.

### 3. Test principle

The wells of the microtiter strips are coated with specific antibodies to egg white proteins. By adding standards and samples to the wells, egg white protein present will bind to the specific antibodies. In a washing step components not bound are removed. Then antibody conjugated to peroxidase is added. This antibody conjugate is bound to the Ab-Ag-complex. An antibody-antigen-antibody (sandwich) complex is formed. Any unbound conjugate is then removed in a washing step. The detection of egg white protein takes place by adding Substrate/Chromogen. The enzyme conjugate converts the chromogen into a blue product. The addition of the stop solution leads to a color change from blue to yellow. The measurement is made photometrically at 450 nm. The absorbance is proportional to the egg white proteins concentration of the sample. The result is expressed in mg/kg whole egg powder.

## 4. Reagents provided

Each kit contains sufficient materials for 48 measurements (including standard analyses).

Component	Cap color	Format		Volume
Microtiter plate	-	Ready to use		48 wells
Allergen Extraction buffer	Green	<b>Concentrate</b>	<b>10x</b>	100 ml
Standard 1*	Transparent	Ready to use	0 mg/kg	1.3 ml
Standard 2*	Transparent	Ready to use	0.5 mg/kg	1.3 ml
Standard 3*	Transparent	Ready to use	1.5 mg/kg	1.3 ml
Standard 4*	Transparent	Ready to use	4.5 mg/kg	1.3 ml
Standard 5*	Transparent	Ready to use	13.5 mg/kg	1.3 ml
Wash buffer	Brown	<b>Concentrate</b>	<b>10x</b>	100 ml
Conjugate	Red	<b>Concentrate</b>	<b>11x</b>	0.7 ml
Substrate/Chromogen Red Chromogen Pro	Brown	Ready to use		10 ml
Stop solution	Yellow	Ready to use		14 ml

- \*) The **dilution factor of 20** for the sample has already been taken into account. Therefore, the whole egg powder concentration of samples can directly be read from the standard curve.

## 5. Materials required but not provided

### 5.1. Equipment:

- microtiter plate spectrophotometer (450 nm)
- centrifuge + centrifugal vials
- shaker
- water bath
- laboratory mincer/grinder, mortar, ultra-turrax or homogenizer
- graduated pipettes
- variable 20 - 200 µl and 200 - 1000 µl -micropipettes

### 5.2. Reagents:

- distilled or deionized water

## 6. Warnings and precautions for the users

This test should only be carried out by trained laboratory employees. The instruction for use must be strictly followed.

This kit may contain hazardous substances. For hazard notes on the contained substances please refer to the appropriate material safety data sheets (MSDS) for this product, available online at [www.r-biopharm.com](http://www.r-biopharm.com).

## 7. Storage instructions

Store the kit at 2 - 8 °C (35 - 46 °F). **Do not freeze any test kit components.**

Return any unused microwells to their original foil bag, reseal them together with the desiccant provided and further store at 2 - 8 °C (35 - 46 °F).

The reddish substrate/chromogen is light sensitive, therefore, avoid exposure to direct light.

No quality guarantee is accepted after the expiration date on the kit label.

Do not interchange individual reagents between kits of different lot numbers.

## 8. Indication of instability or deterioration of reagents

- any bluish coloration of the red stained substrate-/chromogen prior to test implementation
- a value of less than 0.8 absorbance units ( $A_{450\text{ nm}} < 0.8$ ) for standard 5

## 9. Preparation of Samples

Working devices such as a mill, glass vials or spatulas must be cleaned before and after each sample preparation to remove any remains of egg and to avoid contamination.

The **Allergen Extraction buffer** is provided as a 10fold concentrate. Before dilution of the buffer concentrate dissolve any crystals in a water bath at 37 °C (99°F) completely and mix well. After that dilute the heated buffer concentrate 1:10 (1+9) with distilled water before use (i.e. 100 ml buffer concentrate + 900 ml distilled water). The diluted buffer is stable at 2 - 8 °C (35 - 46 °F) for approx. twelve weeks.

### 9.1. Sample preparation

- grind 5 g of the sample carefully and mix thoroughly
- weigh 1 g of the sample and add 20 ml diluted Allergen Extraction buffer (the extraction buffer should already have a temperature of approx. 60 °C (140 °F))
- or add 19 ml of the diluted Allergen Extraction buffer to 1 ml liquid sample (the extraction buffer should already have a temperature of approx. 60 °C (140 °F))
- mix intensively and incubate for 10 min at 60 °C (140 °F), afterwards cool down
- centrifuge: 10 min, at least 2500 g, if possible at 4 °C (39 °F) and/or filter the extract (alternatively 2 ml of the extract can be centrifuged at high speed for 10 min in reaction caps by using a microcentrifuge)
- use 100 µl per well in the assay

### **Remark:**

**For fenugreek, cloves, mustard or celery containing samples, add 1 g casein to 1 g or 1 ml sample to avoid disturbing effects. When analyzing samples with unknown composition, casein may be added as a precaution.**

**The sample extracts can be stored at 2 - 8 °C (35 - 46 °F) for 1 day.**

## **10. Test implementation**

### 10.1. Test preparation

Bring all reagents to room temperature (20 - 25 °C / 68 - 77 °F) before use.

The **conjugate** (bottle with red cap) is provided as an 11fold concentrate. Since the diluted conjugate has a limited stability, only the amount which actually is needed should be reconstituted. Before diluting, the conjugate concentrate should be shaken carefully. For reconstitution, the conjugate concentrate is diluted 1:11 (1+10) in distilled water (e. g. 200 µl concentrate + 2 ml distilled water, sufficient for 2 microtiter strips).

The **wash buffer** is provided as a 10fold concentrate. Before dilution of the buffer concentrate dissolve any crystals in a water bath at 37 °C (99 °F) completely and mix well. After that dilute the heated buffer concentrate 1:10 (1+9) with distilled water before use (i.e. 100 ml buffer concentrate + 900 ml distilled water). The diluted wash buffer is stable at 2 - 8 °C (35 - 46 °F) for approx. four weeks.

### 10.2. Test procedure

Carefully follow the recommended washing procedure. Do not allow microwells to dry between working steps.

Do not use more than 3 strips (24 wells) at a time. In the case of more than three strips, a second uncoated plate (e.g. low binding from Greiner bio-one Cat.-No. 655101) should be used as a pre-plate to avoid a time shift over the microtiter plate. All standards and samples are pipetted into the uncoated plate (at least 150 µl per well) and then quickly transferred to the coated microtiter plate with an 8-channel pipette.

1. Insert a sufficient number of wells into the microwell holder for all standards and samples to be run. Record standard and sample positions.
2. Add 100 µl of each standard or prepared sample in duplicates<sup>0</sup> to the wells, mix gently by shaking the plate manually and incubate for 10 min at room temperature (20 - 25 °C / 68 - 77 °F).
3. Pour the liquid out of the wells and tap the microwell plate upside down vigorously (three times in a row) against absorbent paper to ensure complete removal of liquid from the wells. Fill all the wells with 250 µl wash buffer (see 10.1.) and pour out the liquid again. Repeat two more times.
4. Add 100 µl of the diluted conjugate (see 10.1.) to each well, mix gently by shaking the plate manually and incubate for 10 min at room temperature (20 - 25 °C / 68 - 77 °F).
5. Pour the liquid out of the wells and tap the microwell plate upside down vigorously (three times in a row) against absorbent paper to ensure complete removal of liquid from the wells. Fill all the wells with 250 µl wash buffer (see 10.1.) and pour out the liquid again. Repeat two more times.
6. Add 100 µl of substrate/chromogen to each well. Mix gently by shaking the plate manually and incubate for 10 min at room temperature (20 - 25 °C / 68 - 77 °F) in the dark.
7. Add 100 µl of the stop solution to each well. Mix gently by shaking the plate manually and measure the absorbance at 450 nm. Read within 10 min after addition of stop solution.

## 11. Results

A special software, the RIDA<sup>®</sup>SOFT Win / RIDA<sup>®</sup>SOFT Win.net (Art. Nr. Z9996), is available for evaluation of the RIDASCREEN<sup>®</sup> enzyme immunoassays. The calculation should be done by use of a cubic spline function. The course of the standard curve is shown in the Quality Assurance Certificate enclosed in the test kit.

In comparison with the certificate, higher values of the absorbance ( $A_{450\text{ nm}}$ ) for the calibration curve, especially for the zero standard, may be a result of insufficient washing or egg contamination.

A further dilution and new detection of the samples is recommended for absorbance values ( $A_{450\text{ nm}}$ ) > standard 5.

### **Please note:**

When working according to the described sample preparation, the dilution factor is 20. The allergen concentration can be read directly from the standard curve (the sample dilution factor of 20 is already taken into account for the standard concentrations (see 4. \*)).

Sample dilutions of more than 1:20 must be considered for the calculation of the whole egg powder concentration.

The test is calibrated against NIST Reference Material 8445 (whole egg powder). The result is expressed as mg/kg whole egg powder. The reference material 8445 contains 49% +/- 1% total protein. If the result is multiplied by 0.49, then the result can be expressed as mg/kg (ppm) total protein. If the result is multiplied by 0.263 then mg/kg (ppm) egg white protein are obtained.

### **Calculation example**

From the standard curve 10 mg/kg whole egg powder are read. If multiplied by 0.49 then 4.9 mg/kg (ppm) whole egg protein is obtained. If multiplied by 0.263 then 2.63 mg/kg (ppm) egg white protein is obtained.

### **In general**

Samples tested negative still could contain an allergen contamination below the limit of detection of the assay, or they might contain other allergen components like lipids for example.

Due to the multitude of food types, matrix effects cannot be excluded. In processed food (e.g. heat treatment, dehydration, etc.), proteins may be altered or fragmented, this may have an impact on the recovery/cross reactivity.

Allergen containing samples that have been heat treated show a reduced recovery because the proteins denature and are no longer recognized by the antibody. The reduction in recovery depends strongly on the temperature and the duration of the heat treatment. If samples are heat treated at high temperature the recovery can be significantly reduced.

For evaluation of the cross reactivity only one exemplary sample was analyzed, other samples may show a different result. All cross reactivities and exemplary analyzed matrices are described in the Validation report

## **Recommendation**

In order to ensure a high analytical performance we recommend:

- to analyze each sample material should in duplicates
- to use also allergen-free and allergen- containing (spiked) samples as test controls
- due to the multitude of food types, matrix effects cannot be excluded. To ensure an accurate result spike experiments are recommended
- to perform SureFood<sup>®</sup> PCR for confirmation of the result
- For details using the ChemWell<sup>®</sup> or GEMINI automation please contact [sales@r-biopharm.de](mailto:sales@r-biopharm.de).

For further product information and application notes, please contact [sales@r-biopharm.de](mailto:sales@r-biopharm.de).

**The product information folder with further information is available on request from your local distributor or R-Biopharm AG.**

## **Further application notes**

- Sample preparation for wine
- RIDASCREEN<sup>®</sup>FAST Ei/Egg Protein (Art. No. R6402) used with automates GEMINI and ChemWell<sup>®</sup> 2910 Automate
- RIDASCREEN<sup>®</sup>FAST – Swabbing Allergen Swabbing method for the qualitative analysis of allergens in the production line or for laboratory equipment

The data corresponds to our present state of technology and provides information on our products and their uses. R-Biopharm makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. Defective products will be replaced. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. R-Biopharm shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

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are registered trademarks of R-Biopharm AG