# **Cocktail ECO**

**Art. No. R7080**

In vitro Test

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Storage at 2 - 8 °C

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# Cocktail ECO

## 1. Intended use

The Cocktail ECO can be used for the sample preparation of raw materials, heat treated and processed food.

The Cocktail ECO is used in combination with the following test systems:

- RIDASCREEN® Gliadin (R7001)
- RIDASCREEN®FAST Gliadin (R7002)
- RIDASCREEN®FAST Gliadin sensitive (R7051)
- RIDA®QUICK Gliadin (R7003)
- RIDA®QUICK Gliadin (single packaged) (R7004)

The sample preparation using the **Cocktail (patented)** (R7006/R7016) is the official R5-Mendez method according to the Codex Alimentarius and the AOAC. The faster sample preparation using the environmental-friendly **Cocktail ECO** (R7080) is convenient for the screening of samples. The Cocktail ECO has an extraction efficiency of approx. 70 - 110% compared to Cocktail (patented).

Sample preparation: homogenization and extraction

Time requirement: Cocktail ECO (for 10 samples) ..... approx. 35 min.

## 2. Reagents provided

Component	Cap color	Format	Volume
Buffer ECO	green	Ready to use	2 x 115 ml
Additive ECO	green	Ready to use	6 g

## 3. Warnings and precautions for the users

This test should only be carried out by trained laboratory employees. The instructions for use must be strictly followed.

The Cocktail Eco contains no hazardous substances. It is not necessary to work under a chemical hood. For hazard notes on the contained substances please refer to the appropriate material safety data sheets (MSDS) for this product, available online at [www.r-biopharm.com](http://www.r-biopharm.com).

## 4. Materials required but not provided

### 4.1. Equipment

- centrifuge, centrifugal vials (e.g. Brand 10742512)
- shaker
- laboratory mincer / grinder, pestle and mortar, ultra-turrax or homogenizer
- water bath (50 °C / 122 °F)

### 4.2. Reagents

- distilled or deionized water
- gluten-free **skim milk powder** (food quality)
- **Ethanol solution (80 %)**: i.e. add 120 ml ethanol p.a. to 30 ml distilled water and shake well

## 5. Storage instructions

Store the kit at 2 - 8 °C (35 - 46 °F) (do not freeze).

## 6. Preparation of Samples

### 6.1. Preliminary comments

Airborne cereal dust and dirty laboratory equipment lead to gluten contamination of the assay. Therefore, before starting and during the assay wear gloves.

- clean surfaces, glass vials, mincers and other equipment with 40 % ethanol or 2-propanol
- carry out the sample preparation in a room isolated from the ELISA procedure
- check for gluten contamination of reagents and equipment with the test strips RIDA<sup>®</sup>QUICK Gliadin (R7003/R7004/R7005)

### 6.2. Mixing of the Cocktail ECO

To produce the Cocktail ECO it is necessary to add the Additive ECO to the Buffer ECO. Only mix the quantity of Cocktail ECO that is actually needed because the stability of the Cocktail ECO is only one day. The volume needed depends on the weigh in. The minimal weigh in is 0.25 g. **Due to the fact that the gluten contamination is often very heterogeneous a sample weigh in of 1 g is recommended.** For a sample weigh in of 1 g use 10 ml Cocktail ECO and 30 ml 80 % Ethanol for the following sample preparation.

Number of samples and weigh in	Production Cocktail ECO		Needed volume per sample	
	Buffer ECO	Additive ECO	Cocktail ECO	80% Ethanol
1 sample à 0.25 g	3 ml	60 mg	2.5 ml	7.5
5 samples à 0.25 g	15 ml	300 mg	2.5 ml	7.5
10 samples à 0.25 g	27 ml	540 mg	2.5 ml	7.5
20 samples à 0.25 g	56 ml	1,08 g	2.5 ml	7.5
1 sample à 1 g	12 ml	240 mg	10 ml	30
10 samples à 1 g	120 ml	2,4 g	10 ml	30

### 6.3. Preparations

Homogenize well a sufficient amount (at least 5 g or 5 ml) of sample (grind it thoroughly to powder and mix well or mix well the solution respectively).

- **liquid food samples:** use 0.25 ml of the homogenized sample and add 2.5 ml of the Cocktail ECO, close the vial and mix well
- **other food samples (e.g. soy and quinoa containing samples):** weigh 0.25 g of the homogenized sample and add 2.5 ml of the Cocktail ECO, close the vial and mix well
- **tannin and polyphenol containing food samples (e.g. chocolate, coffee, cocoa, chestnut flour, buckwheat, millet and spices):** weigh 0.25 g of the homogenized sample, add 0.25 g of skimmed milk powder and add 2.5 ml of the Cocktail ECO, close the vial and mix well
- **meat and sausages:** in these matrices gliadin may not be distributed evenly, therefore, homogenized 50 g sample: weigh 0.25 g of the homogenized sample and add 2.5 ml of the Cocktail ECO, close the vial and mix well
- **oat samples:** gliadin may not be distributed evenly, furthermore the samples are difficult to homogenize. Therefore, homogenize 200 g, then carry out the extraction with at least the fourfold amount of reagents: weigh 1 g of the homogenized sample and add 10 ml of the Cocktail ECO, close the vial and mix well

**Please further extract all samples as described in the following:**

- incubate for 10 min at 50 °C (122 °F)
- let the sample cool down and then mix it with 7.5 ml 80 % ethanol (for 1 g sample weigh in: 30 ml 80 % ethanol)
- close the vial and shake for 10 min up side down or by a rotator at room temperature (20 - 25 °C / 68 - 77 °F)

- centrifuge: 5 min, at least 2500 g, at room temperature (20 - 25 °C / 68 - 77 °F) and/or filter the extract (alternatively 2 ml of the extract can be centrifuged with high speed for 10 min in reaction caps by using a microcentrifuge)
- transfer the supernatant in a screw top vial

**Remark:**

The supernatant obtained after the centrifugation or the filtrate can be stored in a tightly closed vial in the dark at room temperature (20 - 25 °C / 68 - 77 °F) up to two weeks.

When using the **ELISA RIDASCREEN® Gliadin (R7001)**, **RIDASCREEN®FAST Gliadin (R7002)** und **RIDASCREEN®FAST Gliadin Sensitive (R7051)** (see also the respective test kit insert):

- dilute the sample 1:12.5 (1+11.5 / 80 µl + 920 µl) with buffer: the final dilution factor is 500 (Attention: for R7001 and R7002 the sample dilution buffer needs to be diluted 1:5 before use, see respective test kit insert)
- use immediately 100 µl per well in the assay.

When using the **lateral flow test strip RIDA®QUICK Gliadin** (see also respective test kit insert):

- 500 µl of the sample dilution buffer from the test kit are transferred to a test tube
- 50 µl of the extracted sample are added.

**Recommendation:**

In order to ensure a high analytical performance:

- Use also gluten free and gluten containing (spiked) samples as test controls
- Carry out spiking experiments for an accurate and correct procedure
- Strong acidic or alkaline samples have to be neutralized prior to analysis
- Confirm results with PCR (e.g. SureFood® PCR Allergen Gluten)
- In the production of foods such as beer or sourdough, proteins are fragmented. In sandwich ELISAs protein fragments lead to a reduced recovery, such samples should be analyzed with a competitive ELISA test systems like the RIDASCREEN® Gliadin competitive (R7021)

Further product information and application notes, please contact [sales@r-biopharm.de](mailto:sales@r-biopharm.de).

The data corresponds to our present state of technology and provides information on our products and their uses. R-Biopharm makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. Defective products will be replaced. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. R-Biopharm shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

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are registered trademarks of R-Biopharm AG
### **R-Biopharm AG**