

PDS No. 6901xx / 6901xx-TRI	PRODUCT DATA SHEET					Page 1 of 2
Revision 11	Cell Culture Flask, CELLSTAR®, 50 ml					 greiner bio-one
	Greiner Item-No. 6901xx / 6901xx-TRI					
Valid for Item-No.:	690160 (sterile)	690170 (sterile)	690175 (sterile)	690190 (sterile)	690195 (sterile)	
	690160-TRI (sterile)		690175-TRI (sterile)	690190-TRI (sterile)	690195-TRI (sterile)	

1.	Description / Specification	
1.1	Description	Cell Culture Flask, 50 ml, canted neck, printed and engraved graduation (scale 5-40 ml), sterile 690160: physical surface treatment, standard screw cap* 690160-TRI: physical surface treatment, standard screw cap*, triple packed 690170: physical surface treatment, measuring grid, standard screw cap* 690175: physical surface treatment, filter cap 690175-TRI: physical surface treatment, filter cap, triple packed 690190: hydrophobic surface (without treatment) , standard screw cap* 690190-TRI: hydrophobic surface (without treatment) , standard screw cap*, triple packed 690195: hydrophobic surface (without treatment), filter cap 690195-TRI: hydrophobic surface (without treatment), filter cap, triple packed * with ventilation position
1.2	Dimensions	Flask with cap: see Customer Drawing 690175, -175-TRI, -195, -195-TRI: pore size of filter membrane: 0.2 µm
1.3	Volume	Max. volume: 50 ml 690160, -160-TRI, -170, -175, -175-TRI: - working volume: 5-10 ml - growth area: 25 cm²
1.4	Material / Resin	Flask: PS (Polystyrene), free of heavy metal Cap: HDPE (High Density Polyethylene), free of heavy metal Filter: PET (Polyethylene Terephthalate), free of heavy metal
1.5	Colour	Flask: clear; print: graduation white, engraved graduation: clear 690170: measuring grid: clear 690160, -160-TRI, -170, -175, -175-TRI: Cap: red 690190, -190-TRI, -195, -195-TRI: Cap: white 690175, -175-TRI, -195, -195-TRI: Filter: white
1.6	Sterilization	690160, -170, -175, -190, -195: SAL 10⁻³ 690160-TRI, -175-TRI, -190-TRI, -195-TRI: SAL 10⁻⁶
1.7	Quality Control	- Raw Material-Control: physical testing - Product-Control: testing of attributive and variable characteristics in accordance with the valid specification
1.8	Other Information	- For single use only - 690160, -160-TRI, -170, -175, -175-TRI: Expiry date and Lot-No. printed on bottom of flask 690190, -190-TRI, -195, -195-TRI: "Suspension" and Lot-No. printed on bottom of flask

2.	Features	
2.1	Basic features	Free of detectable DNase/RNase, human DNA and pyrogens. Contents non-cytotoxic.
2.2	Temperature range	-20°C to +60 °C
2.3	Autoclavability	No
2.4	Centrifugation, max. RCF	N/A
2.5	Chemical Resistance	See homepage: https://www.gbo.com/en_INT/know-how-services/download-center.html
2.6	Shelf life	690160, -160-TRI, -170, -175, -175-TRI: 4 years after month of production 690190, -190-TRI, -195, -195-TRI: 5 years after month of production

Data Sheet subject to change without notice!

Prior Issue	Drawn	Approved	Released	CONFIDENTIAL: Information contained in this document or drawing is confidential and proprietary to Greiner Bio-One GmbH. This document may not be reproduced for any reason without written permission from Greiner Bio-One GmbH. All rights of design, invention, and copyright are reserved.
Revision 10	Date 30 November 2016	Date 20 June 2017	Date 21 June 2017	
Date 18.08.2016	Name S. Kaelberer	Name Dr. L. Breth	Name A. Schulz	

DISCLAIMER: The description of a certain product can only be considered as a guidance, because its performance ultimately depends on what the product is used for. Very often performance studies are indispensable.

PDS No. 6901xx / 6901xx-TRI	PRODUCT DATA SHEET				Page 2 of 2
Revision 11	Cell Culture Flask, CELLSTAR®, 50 ml				 greiner bio-one
	Greiner Item-No. 6901xx / 6901xx-TRI				
Valid for Item-No.:	690160 (sterile)	690170 (sterile)	690175 (sterile)	690190 (sterile)	690195 (sterile)
	690160-TRI (sterile)		690175-TRI (sterile)	690190-TRI (sterile)	690195-TRI (sterile)

3.	Packaging	690160, -170, -175, -190, -195	690160-TRI, -175-TRI, -190-TRI, -195-TRI
3.1	Pieces / Bag	10	10 (triple packed)
3.2	Pieces / Box	200	130
3.3	Lot-No.	E YY MM XXX (manufacturing facility, year, month, consecutive SAP-No.)	
3.4	Other Information	Certificate of Quality	

4.	Other Information
4.1	Laboratory use only

Data Sheet subject to change without notice!

Prior Issue	Drawn	Approved	Released	CONFIDENTIAL: Information contained in this document or drawing is confidential and proprietary to Greiner Bio-One GmbH. This document may not be reproduced for any reason without written permission from Greiner Bio-One GmbH. All rights of design, invention, and copyright are reserved.
Revision 10	Date 30 November 2016	Date 20 June 2017	Date 21 June 2017	
Date 18.08.2016	Name S. Kaelberer	Name Dr. L. Breth	Name A. Schulz	

DISCLAIMER: The description of a certain product can only be considered as a guidance, because its performance ultimately depends on what the product is used for. Very often performance studies are indispensable.