


PDS No. 6041xx / 604160-TRI	<b>PRODUCT DATA SHEET</b>				Page 1 of 1
Revision 10	Pipettes, CELLSTAR®, 1 ml				 greiner bio-one
	Greiner Item-No. 6041xx / 604160-TRI				
Valid for Item-No.:	<b>604107</b> (sterile)	<b>604160</b> (sterile)	<b>604181</b> (sterile)	<b>604160-TRI</b> (sterile)	

1.	Description / Specification	
1.1	Description	PS Pipette, graduation (scale 1/100 ml), lateral yellow stripe, synthetic plug, 1 ml, sterile
1.2	Dimensions	See Customer Drawing
1.3	Volume	1 ml; volumetric accuracy +/- 2 %
1.4	Material / Resin	Pipette: PS (Polystyrene), free of heavy metal Plug: Bonded Polyester, free of heavy metal
1.5	Colour	Pipette: clear Plug: white
1.6	Sterilization	604107, -160, -181: SAL 10 <sup>-3</sup> 604160,-TRI: SAL 10 <sup>-6</sup>
1.7	Quality Control	Product-Control: testing of attributive and variable characteristics in accordance with the valid specification
1.8	Other Information	For single use only

2.	Features	
2.1	Basic features	Free of detectable DNase/RNase, human DNA and pyrogens. Contents non-cytotoxic.
2.2	Temperature range	0°C to +40°C
2.3	Autoclavability	No
2.4	Centrifugation, max. RCF	N/A
2.5	Chemical Resistance	See homepage: <a href="https://www.gbo.com/en_INT/know-how-services/download-center.html">https://www.gbo.com/en_INT/know-how-services/download-center.html</a>
2.6	Shelf life	5 years after month of production
2.7	Other Information	-

3.	Packaging	604107	604160	604181	604160-TRI
3.1	Pieces / Bag	25 Bulk packs (plastic)	1 / 100 "peel off"- packaging (plastic/plastic)	1 / 100 "peel off"- packaging (paper/plastic)	1 "peel off"- packaging (plastic/plastic)  10 "Triple packed"
3.2	Pieces / Box	1000	1000	1000	100
3.3	Lot-No.	F YY MM XXX (manufacturing facility, year, month, consecutive SAP-No.)	F YY MM XXX (manufacturing facility, year, month, consecutive SAP-No.)	F YY MM XXX (manufacturing facility, year, month, consecutive SAP-No.)	E YY MM XXX (manufacturing facility, year, month, consecutive SAP- No.)
3.4	Other Information	Certificate of Quality			

4.	Other Information
	-

Data Sheet subject to change without notice!

Prior Issue	Drawn	Approved	Released	CONFIDENTIAL: Information contained in this document or drawing is confidential and proprietary to Greiner Bio-One GmbH. This document may not be reproduced for any reason without written permission from Greiner Bio-One GmbH. All rights of design, invention, and copyright are reserved.
Revision 09	Date 12 April 2019	Date 15 April 2019	Date 15 April 2019	
Date 21.11.2016	Name S. Kaelberer	Name Dr. T. Schreiber	Name A. Illig	

**DISCLAIMER:** The description of a certain product can only be considered as a guidance, because its performance ultimately depends on what the product is used for. Very often performance studies are indispensable.