


PDS No. 6071xx / 607160-TRI	<b>PRODUCT DATA SHEET</b>				Page 1 of 1
Revision 09	Pipettes, CELLSTAR®, 10 ml				 greiner bio-one
	Greiner Item-No. 6071xx / 607160-TRI				
Valid for Item-No.:	<b>607107</b> (sterile)	<b>607160</b> (sterile)	<b>607180</b> (sterile)	<b>607160-TRI</b> (sterile)	

1.	Description / Specification	
1.1	Description	PS Pipette, graduation (scale 1/10 ml), lateral orange stripe, with synthetic plug, 10 ml, sterile
1.2	Dimensions	See Customer Drawing
1.3	Volume	10 ml; volumetric accuracy $\pm 2\%$
1.4	Material / Resin	Pipette and nozzle: PS (Polystyrene), free of heavy metal Plug: Bonded Polyester, free of heavy metal
1.5	Colour	Pipette: clear Plug: white
1.6	Sterilization	607107, -160, -180: SAL $10^{-3}$ 607160-TRI: SAL $10^{-6}$
1.7	Quality Control	- <u>Raw Material-Control</u> : physical testing - <u>Product-Control</u> : testing of attributive and variable characteristics in accordance with the valid specification
1.8	Other Information	For single use only

2.	Features	
2.1	Basic features	Free of detectable DNase/RNase, human DNA and pyrogens. Contents non-cytotoxic.
2.2	Temperature range	0°C to +40°C
2.3	Autoclavability	No
2.4	Centrifugation, max. RCF	N/A
2.5	Chemical Resistance	See homepage: <a href="https://www.gbo.com/en_INT/know-how-services/download-center.html">https://www.gbo.com/en_INT/know-how-services/download-center.html</a>
2.6	Shelf life	5 years after month of production
2.7	Other Information	-

3.	Packaging	607107	607160	607180	607160-TRI
3.1	Pieces / Bag	25 Reclosable bulk packs (plastic)	1 "peel off"-packaging (plastic/plastic)	1 "peel off"-packaging (paper/plastic)	1 "peel off"-packaging (plastic/plastic)  10 "Triple packed"
3.2	Pieces / Box	500	200	200	100
3.3	Lot-No.	F YY MM XXX (manufacturing facility, year, month, consecutive SAP-No.)	F YY MM XXX (manufacturing facility, year, month, consecutive SAP-No.)	F YY MM XXX (manufacturing facility, year, month, consecutive SAP-No.)	E YY MM XXX (manufacturing facility, year, month, consecutive SAP-No.)
3.4	Other Information	Certificate of Quality			

4.	Other Information
4.1	-

Data Sheet subject to change without notice!

Prior Issue	Drawn	Approved	Released	CONFIDENTIAL: Information contained in this document or drawing is confidential and proprietary to Greiner Bio-One GmbH. This document may not be reproduced for any reason without written permission from Greiner Bio-One GmbH. All rights of design, invention, and copyright are reserved.
Revision 08	Date 30 January 2019	Date 31 January 2019	Date 31 January 2019	
Date 01.10.2015	Name S. Kaelberer	Name Dr. J. Bischoff	Name A. Illig	

**DISCLAIMER:** The description of a certain product can only be considered as a guidance, because its performance ultimately depends on what the product is used for. Very often performance studies are indispensable.